



Conformance Statement
Version 5.5

IHE Integration Statement

VISUS Health IT GmbH
Gesundheitscampus-Süd 15
D-44801 Bochum

www.visus.com

Document type: Conformance Statement
Product: JiveX
Version: 5.5
Creation date: December 6, 2023

Used symbols:



Note:

This symbol indicates special information for easier product operation or it provides other important information.



Warning:

This warning symbol indicates important safety-related information, like warnings and precautions which cannot be placed on the product itself.

Regulatory:

CE 0482 JiveX is a class IIa medical device in accordance with regulation (EU) 2017/745.

Notice to users in the European Union: Any serious incident that has occurred in relation to the medical device should be reported to the manufacturer and the competent authority.

© 2023 VISUS Health IT GmbH

JiveX® is an internationally registered trademark of VISUS Health IT GmbH. All other product names are trademarks or trade names of their respective owners.

Table of contents

1 Supported IHE Profiles 4

1 Supported IHE Profiles

IHE Integration Statement

Date 2023–10–31

Vendor	Product Name	Version
VISUS Health IT GmbH	JiveX	5.4
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information (ARI)	Image Display	None
	Image Manager/Archive	None
Audit Trail and Node Authentication (ATNA)	Secure Application	None
Basic Patient Private Consents (BPPC)	Content Consumer	None
Consistent Presentation of Images (CPI)	Image Display	None
	Image Manager/Archive	None
Consistent Time (CT)	Time Client	None
Digital Breast Tomosynthesis (DBT)	Image Display	None
	Image Manager/Archive	None
Evidence Documents (ED)	Image Display	None
	Image Manager/Archive	None
Import Reconciliation Workflow (IRWF)	Image Manager/Archive	None
	Importer	None
	Performed Procedure Step Manager	None
Invoke Image Display (IID)	Image Display	None
Key Image Note (KIN)	Image Display	None
	Image Manager/Archive	None
Mammography Image (MAMMO)	Image Display	None
	Image Manager/Archive	None
Patient Administration Management (PAM)	Patient Demographics Consumer	None
	Patient Encounter Consumer	None
	Performed Procedure Step Manager	None
Patient Demographics Query (PDQ)	Patient Demographics Consumer	None
Patient Information Reconciliation (PIR)	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Portable Data for Imaging (PDI)	Image Display	DVD, USB
	Portable Media Creator	Web Content, DVD, USB
	Portable Media Importer	DVD, USB

IHE Integration Statement

Date 2023–10–31

Scheduled Workflow (SWF)	Image Display	None
	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Scheduled Workflow.b (SWF.b)	Image Display	None
	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Sharing Value Sets (SVS)	Value Set Repository	None
Cross-Enterprise Document Reliable Interchange (XDR)	Document Recipient	None
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer	None
	Document Repository	None
	Integrated Document Source/Repository	Document Replacement, Document Addendum
	Document Source	Document Replacement, Document Addendum
Cross-enterprise Document Sharing for Imaging.b (XDS-I.b)	Imaging Document Consumer	None
Cross-Enterprise User Assertion (XUA)	X-Service User	None
Internet address for vendor's IHE information		
www.visus.com/en/download/conformance-statements.html		
Links to Standards Conformance Statements for the implementation		
DICOM/HL7	www.visus.com/en/download/conformance-statements.html	
Links to general information on IHE		
North America http://ihe.net/IHE_Worldwide/	Europe http://ihe.net/IHE_Worldwide/	Asia/Pacific http://ihe.net/IHE_Worldwide/