

IHE Integration Statement

Conformance Statement


Version: 5.1.2



VISUS Health IT GmbH
Gesundheitscampus-Süd 15-17
44801 Bochum, Germany
www.visus.com

Creation date: 09.07.2019

Regulatory:



JiveX is a class IIa medical device within the scope of directive 93/42/EEC.



JiveX has Section 510(k) clearance under the following number: K181964.

© 2019 VISUS Health IT GmbH

JiveX® is an internationally registered trademark. All rights reserved by VISUS Health IT GmbH Germany. Licensees and distribution partners receive a limited usage right.

Table of Content

1 Supported IHE Profiles 4

1 Supported IHE Profiles

IHE Integration Statement		Date
		2019-07-09
Vendor	Product Name	Version
VISUS Health IT GmbH	JiveX	5.1.2
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Image Display	None
	Image Manager/Archive	None
Audit Trail and Node Authentication	Secure Application	None
BPPC (Basic Patient Private Consents)	Content Consumer	None
Consistent Presentation of Images	Image Display	None
	Image Manager/Archive	None
Consistent Time	Time Client	None
Digital Breast Tomosynthesis	Image Display	None
	Image Manager/Archive	None
Evidence Documents	Image Display	None
	Image Manager/Archive	None
Import Reconciliation Workflow	Image Manager/Archive	None
	Importer	None
	Performed Procedure Step Manager	None
Invoke Image Display	Image Display	None
Key Image Note	Image Display	None
	Image Manager/Archive	None
Mammography Image	Image Display	None
	Image Manager/Archive	None
Patient Administration Management	Patient Demographics Consumer	None
	Patient Encounter Consumer	None
	Performed Procedure Step Manager	None
Patient Demographics Query	Patient Demographics Consumer	None
(continued on next page)		

IHE Integration Statement		Date
		2019-07-09
Vendor	Product Name	Version
VISUS Health IT GmbH	JiveX	5.1.2
(continued from previous page)		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Patient Information Reconciliation	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Portable Data for Imaging	Image Display	DVD, USB
	Portable Media Creator	Web Content, DVD, USB
	Portable Media Importer	DVD, USB
Scheduled Workflow	Image Display	None
	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Scheduled Workflow.b	Image Display	None
	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Sharing Value Sets	Value Set Repository	None
XDS (Cross Enterprise Document Sharing)	Document Consumer	None
XDS-I (Cross Enterprise Document Sharing for Images)	Imaging Document Consumer	None
XUA (Cross Enterprise User Assertion)	X-Service User	None
Internet address for vendor's IHE information:		
www.visus.com/en/download/conformance-statements.html		
Links to Standards Conformance Statements for the implementation		
DICOM/HL7	www.visus.com/en/download/conformance-statements.html	
Links to general information on IHE		
North America: http://ihe.net/IHE_Worldwide/	Europe: http://ihe.net/IHE_Worldwide/	Asia/Pacific: http://ihe.net/IHE_Worldwide/